



# FTS-400 Floor-to-Stand For Home and Care Centers

**The FTS-400 was developed for mobility challenged individuals and those who care for them!**

## Description

The FTS-400 is a lift & transfer tool that enables quick and safe fall recovery while eliminating manual lifting.

It can raise individuals up to 400 lbs. (181 kg) to a seated or standing position, and allows for easy transfers to and from a wheelchair, bed, commode, exercise equipment or even the floor.

The standing lift function allows users to be safely raised to access a high bed or an overhead cabinet - a step stool without the step!

Can be used independently or with the assistance of a caregiver.




## Assisted Operation



## Self Operation




**Standing  
Function**



# Summary of Features

- Beveled seat-front for easy slide on transfer surface
- Wired remote for convenient operation
- Swinging/Removable Rise-Assist Handles for easy secure transfer
- Dual safety belts for added safety
- Powder coated for easy cleaning
- Small footprint for tight spaces
- Lithium-Ion battery powered with automatic charger
- 9' (2.74m) power cord



## Physical Specifications

Safe Working Load .....	400 lbs. (181kg)
Overall Depth .....	26.5 " (673mm)
Maximum Overall Height (raised) .....	70.75" (1746mm)
Minimum Overall Height (lowered) ...	41.75" (1009mm)
External Width incl Rise Handles.....	23" (584mm)
Maximum Seat Height .....	30" (762mm)
Seat Width .....	14" (356mm)
Wheel Diameter.....	5" (127mm)
Unit Weight .....	78 lbs. (35kg)

## Electrical Specifications

Standard 2-prong 110 VAC Power Cord (USA)  
Standard 2-prong 220 VAC Power Cord (International)  
Operating Environment: >+5° C <+40° C (41° to 104°F)  
Lithium-Ion Battery: 24 VDC

## Maintains Privacy, Dignity and Independence



[info@callbeforeyoufall.com](mailto:info@callbeforeyoufall.com)

<https://callbeforeyoufall.com/lifting-chairs/>

800-829-1491